



### Minimally invasive treatment of skin lesions

The CryoPro® is a handheld medical device for cutaneous cryosurgery that uses localized application of liquid nitrogen to treat skin lesions by tissue freezing.

The CryoPro® delivers a continuous flow of liquid nitrogen in a safe, precise, and controlled manner using either open spray or contact probe techniques.

### Simplicity and reliability

The CryoPro® has a long history of reliability, safety, and ease of use and has several advantages compared to traditional excisional surgery:

- Good cosmetic results
- Lower risk of infection
- Minimal wound care
- Low cost

### What's included?

Each CryoPro® device comes with a bent spray and four different open spray apertures (A+B+C+D) suitable for the most common skin lesions. In addition, a stand for the tips and a cleaning adapter are included.

Optionally, contact probes are available in various sizes to better fit different lesion geometries and control the lateral spread of the freeze.

Furthermore, Cortex Technology offers solutions for liquid nitrogen handling and storage.

## Our cool CryoPro®

For liquid nitrogen-based cryosurgery

# 350 ml

CryoPro® MINI  
12-14 hrs. static holding time

# 500 ml

CryoPro® MAXI  
20-24 hrs. static holding time

### CryoPro®

A versatile low-cost treatment of common skin lesions

CryoPro®



**As cool as it gets –**

CryoPro® delivers safety, reliability,  
and precision at -196°C



**Cortex**

Niels Jernes Vej 6B  
9220 Aalborg  
Denmark  
+45 9857 4100  
cortex@cortex.dk  
www.cortex.dk